

## Acceptable Line Cleaners

Bio-Pure Evacuation System Cleaner: *(updated 4/28/2021)*

- pH = 6-8 s.u.

Link: [Bio-Pure | Safco Dental Supply](#)

Bio-Pure®

Sable Industries (Made in USA.) ★★★★★ Write a review



Citrizyme: *(updated 10/4/2021)*

- pH = 7.5 s.u.

Link: [Citrizyme® Concentrated Enzymatic Evacuation System Cleaner - Pascal Company Inc \(pattersondental.com\)](#)



**Citrizyme® Concentrated Enzymatic Evacuation System Cleaner - Pascal Company Inc**

Citrizyme® Concentrated Enzymatic Evacuation System Cleaner is also used as an instrument pre-soak and ultrasonic cleaning solution.

- Premium non-foaming enzymatic cleaner
- Fast-acting dual enzymes remove protein and biofilms quickly and leave no residue
- Fresh citrus scent and color

Evacuation System Cleaner: *(updated 4/28/2021)*

- pH = 7.1 s.u.

Link: [Evacuation System Cleaner \(darbydental.com\)](#)



Item #: 951-2111

**Evacuation System Cleaner**

★ Favorite Add to list

Select your item:

951-2111 Cleaner, 64 oz.

Manufacturer: Darby Dental Supply

Catalog Page: 343

Darby's Evacuation System Cleaner is an enzyme-based cleaner that provides residual over-night debris removal from cuspidor drains and suction systems. It is environmentally friendly, foam free, and pleasantly scented. It works under extreme conditions, including cold water, is pH neutral and safe for both equipment and staff. Darby's Evacuation System Cleaner is also safe in amalgam separators.

---

## Monarch CleanStream Evacuation

System Cleaner: *(updated 4/28/2021)*

- pH = 7-8 s.u.

Link: [Monarch CleanStream | Safco Dental Supply](#)

- The unique formulation, designed for daily cleaning, penetrates and helps remove line buildup while maintaining vacuum performance and suction
- Cleans and deodorizes evacuation lines for all wet and dry vacuum systems
- Non-foaming, non-enzymatic concentrated formula is very economical and lasts up to 5 days if pre-mixed
- Compliant with BMP for dental amalgam waste: Non-oxidizing and a pH of 7-8
- Non-corrosive and biodegradable



---

NeutraVac: *(updated 4/28/2021)*

- pH = 7-8 s.u.

Link: [Neutra Vac | Safco Dental Supply](#)

NeutraVac™

Biotrol (Made in USA.) ★★★★★ Write a review



### Product Details

Dual action, super-concentrated vacuum line cleaner dissolves and removes organic and inorganic debris. Neutral pH balanced formula compatible with amalgam separators. Pleasing mountain-fresh scent. Dilutes ½ ounce to 1 quart.

---

OR-Evac: *(updated 4/7/2021)*

- pH = 7.1 s.u.

Link: [OR-Evac Evacuation System Cleaner Concentrated Liquid 32 oz Ea](#) ([henryschein.com](http://henryschein.com))



### OR-Evac Evacuation System Cleaner Concentrated Liquid 32 oz Ea

1026316 | Henry Schein Inc. - 102-6316

#### DESCRIPTION:

OR-Evac Evacuation System Cleaner Concentrated Liquid Lemon 32 oz Each

#### PACKAGING OPTIONS (1)

Login For Price

EA  
EACH

Add

[MSDS/SDS](#)

---

---

PUREVAC SC:  
(updated 4/7/2021)

- pH ~ 7 s.u.

Link: [System Cleaners - Henry Schein Dental](#)



**PUREVAC SC Cleaner Evacuation System 5 Liter Ea**  
3120099 | Dentsply Sirona Preventive - 21135

**DESCRIPTION:**

PUREVAC SC Cleaner Evacuation System 5 Liter Each

**PACKAGING OPTIONS (1)**

Login For Price

EA  
EACH

 Hazardous materials are non-returnable.

---

Safco Concentrated Enzymatic Evacuation Line: (updated 4/28/2021)

- pH = 7.1 s.u.

Link: [Safco Concentrated Enzymatic Evacuation Line \(safcodental.com\)](#)

Safco Concentrated Enzymatic Evacuation Line

Safco (Made in USA.)  Write a review



Product Details

Part of the Safco KleenLine™ Series

Triple enzyme evacuation line cleaner provides exceptional cleaning and deodorizing for all types of systems.

- Eco-friendly with a neutral pH
- Non-foaming.
- Phosphate-free. Safe for use with amalgam separators.
- Sweet clean fragrance.


---

Safco Concentrated Evacuation Line: (updated 4/28/2021)

- pH = 7.1 s.u.

Link: [Safco Concentrated Evacuation Line Cleaner \(safcodental.com\)](#)

Safco Concentrated Evacuation Line Cleaner  
(previously called KleenLine)

Safco (Made in USA.)  Write a review



Product Details

Part of the Safco KleenLine™ Series

Anti-microbial eco-friendly evacuation line cleaner for everyday use in all types of evacuation systems.

- Reduces biomass buildup.
  - Non-foaming, with neutral pH.
  - Phosphate-free. Safe for use with amalgam separators.
  - Convenient metered dispenser bottle.
  - Sweet clean fragrance.
-

---

Sani-Treet Green: *(updated 4/7/2021)*

- pH = 6-8 s.u.

Link: [Sani-Treet Green | Safco Dental Supply](#)

## Sani-Treet® Green

Enzyme Industries (Made in USA.) [☆☆☆☆☆Write a review](#)



### Product Details

Enzymatic concentrate with new generation surfactants that biodegrade quickly and efficiently. Safe for use with any amalgam separator. Nontoxic, non-caustic. Lemongrass-lavender scent.

Dilute 1 ounce into 1 quart of water; Each gallon bottle makes 128 one-quart (32 oz.) treatments.

---

Sani-Treet Plus: *(updated 4/7/2021)*

- pH = 6-8 s.u.

Link: [Sani-Treet Plus | Safco Dental Supply](#)

## Sani-Treet Plus®

Enzyme Industries (Made in USA.) [☆☆☆☆☆Write a review](#)



### Product Details

Non-foaming enzymatic pre-soak solution. Safe for use in ultrasonic cleaners and as an evacuation system cleaner. Biodegradable. Dilute 1 ounce into 1 quart warm water.

1 gallon makes 128 one-quart (32 oz.) treatments.

---

Sani Vac: *(updated 4/28/2021)*

- pH = 7.8 s.u.

Link: [Sani Vac | Safco Dental Supply](#)

## Sani Vac™

Crosstex (Made in Canada.) [☆☆☆☆☆Write a review](#)



### Product Details

Low-foaming multi-enzyme formula with a pH of 7.8 meets the EPA requirement for use with amalgam separators. Pleasant lemon scent.

---

Vac Attak Green: *(updated 4/28/2021)*

- pH = 6-8 s.u.

Link: [Vac Attak Green | Safco Dental Supply](#)

## Vac Attak™ Green

Premier (Made in USA.) ★★★★★ 1 Review



---

Vacuosol Neutral: *(updated 4/28/2021)*

- pH = 7-8 s.u.

Link: [Vacuosol Neutral | Safco Dental Supply](#)

## Vacuosol™ Neutral

Biotrol (Made in USA.) ★★★★★ Write a review



---

Vacuum Clean Tablets: *(updated 4/28/2021)*

- pH = 7 s.u.

Link: [Vacuum Clean | Safco Dental Supply](#)

## Vacuum Clean

Palmero (Made in USA.) ★★★★★ 1 Review



---

Vacuum Shock Tablets: *(updated 4/28/2021)*

- pH = 7 s.u.

Link: [Vacuum Shock | Safco Dental Supply](#)

## Vacuum Shock

Palmero (Made in USA.) ★★★★★ 1 Review





## Unacceptable Line Cleaners (Bleach cleaners, cleaners outside a pH range of 6-8 s.u.)

All-In-One Enzymatic Cleaner: (updated 4/28/2021)

- pH = 4.2 s.u.

Link: [All-In-One Enzymatic Cleaner \(darbydental.com\)](http://darbydental.com)



Item #: 955-1800

### All-In-One Enzymatic Cleaner

[★ Favorite](#) [Add to list](#)

Select your item:

955-1800 Enzymatic Cleaner, Gallon

Manufacturer: Darby Dental Supply  
Catalog Page: 343

Darby's All-In-One Enzymatic Cleaner is an environmentally friendly solution that contains no harsh chemicals. This non-foaming solution will deodorize the evacuation system from suction tip to vacuum pump including traps, eliminating odor on contact. Darby's All-In-One Enzymatic Cleaner can also be used as an instrument pre-soak and ultrasonic solution, one gallon of cleaner makes 32 gallons of solution.

#### Features and Benefits

- Non-toxic

Precision Cleanse Plus Evacuation System: (updated 4/28/2021)

- pH = 7-9 s.u.

Link: [Precision Cleanse Plus Evacuation System Cleaner Liquid Unit Dose 64 oz Ea - Henry Schein Dental](http://Henry Schein Dental)



### Precision Cleanse Plus Evacuation System Cleaner Liquid Unit Dose 64 Ea

5646008 | Midmark Corporation - 77000570

#### DESCRIPTION:

Precision Cleanse Plus Evacuation System Cleaner Liquid Unit Dose 64 oz Each

#### PACKAGING OPTIONS (1)

Login For Price

EA  
EACH

 - +

Add

CASE QUANTITY: 6

Made in the U.S.A. This product is Made in the U.S.A.

MORE PRODUCT INFO [MSDS/SDS](#)

SlugBuster: (updated 4/7/2021)

- pH = 8-8.5 s.u.

Link: [SlugBuster | Safco Dental Supply](http://Safco Dental Supply)

## SlugBuster™

Ramvac (DentalEZ) (Made in USA.) 1 Review



---

Turbo-Vac: (updated 4/28/2021)

- pH = 1.4 +/- 0.3

## Turbo-Vac®

Kerr TotalCare (Made in USA.) ☆☆☆☆☆ Write a review

Link: [Turbo-Vac | Safco Dental Supply](#)



### Product Details

Neutral pH line flush solution dissolves debris, cleans, and deodorizes. Biodegradable, non-foaming superconcentrate. Dispenser mixes air and Turbo-Vac solution to create a powerful mist that dislodges debris.

Dilution ratio: 1 oz. to 1 quart water; Each gallon bottle makes 128 one-quart (32 oz.) treatments.

---

Vac Attak: (updated 4/8/2021)

- pH = 10 s.u.

Link: [Premier Dental | Vac Attak™ Evacuation System Cleaner \(premierdentalco.com\)](#)

Vac Attak



\$67.25

SKU: 9011101

---

Vacu Blast: (updated 4/28/2021)

- pH = 8.5 s.u.

## Vacu Blast™

Premier (Made in USA.) ☆☆☆☆☆ Write a review

Link: [Vacu Blast | Safco Dental Supply](#)



### Product Details

Powerful enzymatic liquid evacuation system cleaner compatible with all amalgam separators. Neutral pH, non-corrosive, non-foaming. Pleasant eucalyptus scent.

Each 32 oz. bottle makes 64 one-quart (32 oz.) treatments.

---



## No Information on pH

Solmetex PowerScrub: (updated 4/28/2021)

- pH = Neutral

Link: [PowerScrub | Safco Dental Supply](#)

PowerScrub™

Solmetex (Made in USA.) ★★★★★ Write a review



VacuKleen: (updated 4/28/2021)

- pH = No information provided

Link: [VacuKleen E2 | Safco Dental Supply](#)

VacuKleen E<sup>2</sup>™

Kulzer (Made in USA.) ★★★★★ Write a review



### Product Details

Powerful protease enzymes allow for easy dissolving of residue in evacuation lines. Simply add contents of pouch to warm water and aspirate through evacuation hoses. Lemon scent. Unit dose packaging.